

**KL-LED.STZ4/STZ4 (dual head)** is suitable for the requirements of various operation conditions and it is ideal operating lighting equipment for modern operation rooms and clean operation rooms.

### Features:

- a. The removable handle sheath can be disinfected at the high temperature of 135 °C.
- b. The no welding titanium alloy arm, Imported spring arm is lifted and lowered easily and freely.
- c. Imported LED cold light source is adopted , with the ultrathin optical lens.
- d. Digital means are adopted for stepless regulation of LED brightness.
- e. Light field size can be adjusted by the handle.
- f. The service life of LED shadowless lamp is long (50,000h).
- g. Color temp can be adjusted in three steps.
- h. The imported switch power supply is adopted to control the voltage, making the work voltage permanently stable.



### Specification:

- (1) Illumination(lux): 40,000-140,000,satellite: 40,000-140,000
- (2) Color temperature(k): 4500±500
- (3) Input power(VA): 90±10%
- (4) It has 72 LED bulbs one lamp holder
- (5) Average bulb life is 50,000hours
- (6) Ceiling height for installation(cm): 270-320
- (7) Imported LED cold light source is adopted , with the ultrathin optical lens.
- (8) Digital means are adopted for stepless regulation of LED brightness.
- (9) Light field size can be adjusted by the handle.
- (10) The removable handle sheath can be disinfected at the high temperature of 135 °C.
- (11) Color temp can be adjusted in three steps.
- (12) The imported switch power supply is adopted to control the voltage, making the work voltage permanently stable.
- (13) The no welding titanium alloy arm, Imported spring arm is lifted and lowered easily and freely.

