


LIST OF AIMTECH DRESSING INSTRUMENT SET


No.	Deskripsi	Tipe / Kode	Qty
1	Surgical scissors, Fig. 2, straight, 14 cm, sharp/blunt	11-002-14	1
2	Dressing forceps delicate, straight, 14 cm	13-002-14	1
3	Dissecting forceps 1x2, 14.5 cm	13-100-14	1
4	Pean forceps, straight, 16 cm	15-000-16	1
5	Ochsner-Kocher forceps, straight, 16 cm	15-100-16	1
6	Crile-Wood needle holder, 15 cm	17-014-15	1
7	Probe,	10-602-14	1

Dressing Instrument Set




Surgical Scissors (Straight)

11-002-14 14,0 cm, 5 1/2 "




Dressing Forceps (Straight)

13-002-14 14,5 cm, 5 3/4 "




Dissecting Forceps

13-100-14 14,5 cm, 5 3/4 "



Pean Forceps (Straight)

15-000-16 16,0 cm, 6 1/4 "



Ochsner-Kocher Forceps (Straight)

15-100-16 16,0 cm, 6 1/4 "



Crile-Wood Needle Holder

17-014-15 15,0 cm, 6 "



Probe, Ø 2 mm

10-602-14 14,5 cm, 5 3/4 "