

LIST OF AIMTECH SECTIO CAESAREAN INSTRUMENT SET

No.	Deskripsi	Tipe / Kode	Qty
1	Scalpel handle, standard No. 3	10-103-01	1
2	Scalpel handle, standard No. 4	10-104-01	1
3	Metzenbaum scissors, straight, Fig. 1, 18 cm	11-200-18	1
4	Metzenbaum scissors, curved, Fig. 1, 18 cm	11-201-18	1
5	Littauer stitch scissors, straight, 14 cm	11-706-14	2
6	Standard dressing forceps, straight, 20 cm	13-000-20	1
7	Dissecting forceps 1x2, 20 cm	13-100-20	1
8	Kelly forceps, straight, 14 cm	15-052-14	2
9	Kelly forceps, curved, 14 cm	15-053-14	2
10	Ochsner-Kocher forceps, curved, 16 cm	15-101-16	2
11	Allis tissue forceps, 5x6, straight, 15 cm	15-766-15	1
12	Babcock tissue forceps, straight, 16 cm	15-760-16	1
13	Mikulicz forceps, curved, 20.5 cm	15-339-20	6
14	Forster-Ballenger forceps, straight, 25 cm	15-736-25	2
15	Backhaus towel forceps, 13 cm	15-801-13	6
16	Mayo-Hegar needle holder, 14 cm	17-002-14	2
17	Mayo-Hegar needle holder, 18 cm	17-002-18	1
18	Volkman retractor, 4-prongs, sharp, 22 cm	16-001-04	2
19	Middeldorpf retractor, 20x22 mm, 22 cm	16-035-02	2
20	Doyen retractor, 53 x 80 mm, 25 cm	16-038-25	1
21	Balfour abdominal retractor, 20 cm spreading 200 mm, lateral blades 100 x 35 mm, central blade 100 x 60 mm	16-277-20	1
22	Yankauer suction tube, 27 cm with large central bore d 4.8 mm	10-911-01	1
23	Grave vaginal specula Fig. 1, 75 x 20 mm	35-013-01	1
24	Grave vaginal specula Fig. 2, 95 x 35 mm	35-013-02	1
25	Grave vaginal specula Fig. 3, 115 x 35 mm	35-013-03	1
26	Instrument tray, 300 x 175 x 30 mm	51-120-30	1

Sectio Caesarean Instrument Set




Scalpel Handle Std. 3

10-103-01 12,0 cm, 4 ¾ "




Scalpel Handle Std. 4

10-104-01 13,5 cm, 5 ¼ "



Metzenbaum Scissor (Straight)

11-200-18 18,0 cm, 7 "




Metzenbaum Scissor (Curved)

11-201-18 18,0 cm, 7 "




Litauer Stitch Scissor (Straight)

11-706-14 14,0 cm, 5 ½ "




Dressing Forceps

13-000-20 20,0 cm, 8 "




Dissecting Forceps

13-100-20 20,0 cm, 8 "




Kelly Forceps (Straight)

15-052-14 14,0 cm, 5 ½ "



Kelly Forceps (Curved)

15-053-14 14,0 cm, 5 ½ "



Ochner-Kocher Forceps (Curved)

15-101-16 16,0 cm, 6 ¼ "



Allis Tissue Forceps (Straight)


15-766-15 15,0 cm, 6 "



Babcock Tissue Forceps


15-760-16 16,0 cm, 6 ¼ "

Sectio Caesarean Instrument Set




Mikulicz Forceps (Curved)

15-339-20 20,5 cm, 8 "



Förster-Ballenger Forceps (Straight)

15-736-25 25,0 cm, 9 ¾ "



Backhaus Towel Forceps

15-801-13 13,0 cm, 5 ⅞ "




Mayo-Hegar Needle Holder

17-002-14 14,0 cm, 5 ½ "




Mayo-Hegar Needle Holder

17-002-18 18,0 cm, 7 "



Volkman Retractor, 4-Prongs

16-001-04 22,0 cm, 8 ¾ "



Middeldorpf Retractor


16-035-02 22,0 cm, 8 ¾ "



Doyen Retractor

16-038-25 25,0 cm, 10 "

max. spread 200 mm




Central valve,
100 x 60 mm

Lateral valves,
100 x 35 mm

Balfour Abdominal Retractor

16-277-20 20,0 cm, 8 "



4,8 mm
10 mm

Yankauer Suction Tube

10-911-01 27,0 cm, 10 ¼ "




lateral fixation
screw

Grave Vaginal Specula

35-013-01 75 x 20 mm

Sectio Caesarean Instrument Set



lateral fixation
screw

lateral fixation
screw

Grave Vaginal Specula

35-013-02 95 x 35 mm

Grave Vaginal Specula

35-013-03 115 x 35 mm



Instrument Trays

51-120-30 300 x 175 x 30 mm