

LIST OF AIMTECH THORACOTOMY INSTRUMENT SET

No.	Deskripsi	Tipe / Kode	Qty
1	Gross-Maier sponge forceps, curved, 27 cm	15-703-27	4
2	Gross-Maier sponge forceps, straight, 27 cm	15-702-27	3
3	Allis tissue forceps, 5x6, straight, 15 cm	15-766-15	2
4	Backhaus towel forceps, 11 cm	15-801-11	8
5	Backhaus towel forceps, 11 cm, with holder for tubing up to d 15 mm	15-803-11	1
6	Scalpel handle, standard No. 3	10-103-01	1
7	Scalpel handle, standard No. 4	10-104-01	2
8	Scalpel handle No. 4 L	10-104-04	1
9	Scalpel remover	10-199-00	1
10	Standard dressing forceps, straight, 14.5 cm	13-000-14	3
11	Standard dressing forceps, straight, 20 cm	13-000-20	2
12	Dissecting forceps 1x2, 14.5 cm	13-100-14	3
13	Dissecting forceps 1x2, 18 cm	13-100-18	2
14	TC Potts-Smith dressing forceps 15 cm	14-064-15	2
15	Taylor, forceps, straight, 17 cm with scarper	13-058-17	1
16	Metzenbaum scissors, curved, Fig. 1, 18 cm	11-201-18	2
17	Metzenbaum scissors, curved, Fig. 1, 23 cm	11-201-23	1
18	Surgical scissors, Fig. 2, straight, 16 cm, sharp/blunt	11-002-16	2
19	TC Mayo-Stille scissors, curved, 17 cm	12-171-17	1
20	Surgical scissors curved, delicate, Fig.1, 12 cm blunt/blunt	11-011-12	1
21	Pean forceps, delicate, straight, 13 cm	15-010-13	12
22	Pean forceps, delicate, curved, 13 cm	15-011-13	12
23	Halsted-Mosquito forceps, curved, 1x2, 12 cm	15-161-12	4
24	Ochsner-Kocher forceps, straight, 20 cm	15-100-20	2
25	Pean forceps, straight, 20 cm	15-000-20	2
26	Gemini forceps, curved, 23 cm	15-237-23	2
27	Gemini forceps, curved, 25 cm	15-237-25	2
28	Mikulicz forceps, curved, 20 cm	15-343-20	4
29	Kocher retractor, 4-prongs, blunt, 22 cm	16-007-04	2
30	Kocher retractor, 6-prongs, blunt, 22 cm	16-007-06	2
31	Langenbeck retractor, 65x20 mm, curved	16-020-01	2
32	Roux retractor, Fig. 1, 14 cm	16-050-01	1
33	Martin retractor, 105 x 27 mm, 25 cm	16-048-02	2
34	Mikulicz retractor, 120 x 50 mm, 26 cm	16-028-03	2
35	Abdominal spatula, malleable, 330 x 30 mm	16-072-03	1
36	Mayo-Hegar TC-needle holder, 20 cm	18-002-20	3

37	De Bakey TC-needle holder, 23 cm	18-004-23	1
38	Ligature needle, right, blunt, 20 cm	17-731-20	2
39	Konig ligature conductor, 3 mm, 19.5 cm	17-765-19	1
40	Konig ligature conductor, 5 mm, 19.5 cm	17-767-19	1
41	Probe, d 2 mm, 16 cm	10-602-16	1
42	Doyen director, curved, 14 cm	10-613-14	1
43	Gross-Maier sponge forceps, curved, 20 cm	15-703-20	1
44	Poole suction cannula, d 10 mm, 22 cm	10-914-10	1
45	Redon guide needle curved, 12 Charr. with step formation	10-705-12	1
46	Redon guide needle curved, 14 Charr. with step formation	10-705-14	1
47	Round bowls, d 116 mm, H 50 mm, 0.35 l	51-103-11	2
48	Needle case perforated, with compartments 150 x 95 x 13 mm	56-112-02	1

Thoracotomy Instrument Set



**Gross-Maier
Forceps (Curved)**

15-703-27 27,0cm, 10³/₄ "




**Gross-Maier
Forceps (Straight)**

15-702-27 27,0cm, 10³/₄ "



Allis Forceps (Straight)

15-766-15 15,0 cm, 6 "



Backhaus Towel Forceps

15-801-11 11,0 cm, 4³/₈ "



Backhaus Towel Forceps

15-803-11 11,0 cm, 4³/₈ "



Scalpel Handle Std.3

10-103-01 12,0 cm, 4³/₄ "




Scalpel Handle Std.4

10-104-01 13,5 cm, 5¹/₄ "



Scalpel Handle Std.4L


10-104-04 21,0 cm, 8¹/₄ "



Packing unit: 100
pieces


**Scalpel Blade
Remover**

10-199-00




**Dressing Forceps
(Straight)**

13-000-14 14,5 cm, 5³/₄ "



**Dressing Forceps
(Straight)**

13-000-20 20,0 cm, 8 "



Dissecting Forceps

13-100-14 14,5 cm, 5³/₄ "

Thoracotomy Instrument Set



Dissecting Forceps

13-100-18 18,0 cm, 7 "




Potts-Smith Dressing Forceps (Straight)

14-064-15 15,0 cm, 6 "




Taylor Forceps (Straight)

13-058-17 17,0 cm, 6 ¾ "




Metzenbaum Scissors (Curved)

11-201-18 18,0 cm, 7 "



Metzenbaum Scissors (Curved)

11-201-23 23,0 cm, 9 "




Surgical Scissors (Straight)

11-002-16 16,0 cm, 6 ¼ "




Mayo-Stillé Scissors, (Curved)

12-171-17 17,0 cm, 6 ¾ "




Surgical Scissors (Curved)

11-011-12 12,0 cm, 5 ⅛ "



Peen Forceps (Straight)

15-010-13 13,0 cm, 5 ⅛ "




Peen Forceps (Curved)

15-011-13 13,0 cm, 5 ⅛ "



Halsted-Mosquito Forceps (Curved)


15-161-12 12,5 cm, 5 "



Ochsner-Kocher Forceps (Straight)


15-100-20 20,0 cm, 8 "

Thoracotomy Instrument Set




Pean Forceps (Straight)

15-000-20 20,0 cm, 8 "




Gemini Forceps (Curved)

15-237-23 23,0 cm, 9 "




Gemini Forceps (Curved)

15-237-25 25,0 cm, 9 ¾ "




Mikulicz Forceps (Curved)

15-343-20 20,0 cm, 8 "



Kocher Retractor 4-prongs

16-007-04 22,0 cm, 8 ¾ "




Kocher Retractor 6-prongs

16-007-06 22,0 cm, 8 ¾ "



Langenbeck Retractor

16-020-01 21,0 cm, 8 ¼ "



Roux Retractor


16-050-01 14,0 cm, 5 ½ "



105 x 27 mm

Martin Retractor

16-048-02 25,0 cm, 10 "



16-028-03 120 x 50 mm

Mikulicz Retractor


16-028-03 26,0 cm, 10 ½ "



malleable

Abdominal Spatula

16-072-03 330 x 30 mm



Mayo-Hegar Needle Holder


18-002-20 20,0 cm, 8 "

Thoracotomy Instrument Set



De Bakey Needle Holder

18-004-23 23,0 cm, 9 "




Ligature Needle (Right)

17-731-20 20,0 cm, 8 "



König Ligature Conductor

17-765-19 19,0 cm, 7 ½ "



König Ligature Conductor

17-767-19 19,0 cm, 7 ½ "



Probe, Ø 2 mm

10-602-16 16,0 cm, 6 ¼ "




Doyen Director (Curved)

10-613-14 14,5 cm, 5 ¾ "




Gross-Maier Forceps (Curved)

15-703-20 20,0 cm, 8 "




Poole Suction Canula

10-914-10 22,0 cm, 8 ⅝ "



Redon Guide Needle (Curved)

10-705-12 16,0 cm, 6 ¼ "



Redon Guide Needle (Curved)

10-705-14 16,0 cm, 6 ¼ "



Round Bowls

51-103-11 ; Ø 116 x 50 mm; 0,35 l



Needle Case Compartments

56-112-02 150 x 95 x 13 mm