

AKL: 21603914924

## Portable Mono-polar Electrosurgical Unit

### 1 Features

- Has 2 options of output control modes:  
Foot Control Mode and Hand Control Mode
- Mono-polar Cutting and Mono Polar Coagulation
- The output of each watts is manipulated precisely, so that it is applicable to all kinds of micro surgery.



### 2 Output Mode

	Rated output power	Effect
Pure Cut	(1W~100W, Load: 500 $\Omega$ )	The coagulation is effected at the instant of startup, and it's applicable to micro surgery
Blend Cut	(1W~50W, Load: 500 $\Omega$ )	Only with a slight coagulation in the cutting strong
Coagulation	(1W~80W, Load: 500 $\Omega$ )	It is also called strong coagulation or spray coagulation, it coagulates comparatively much faster and the working area is much more special
Soft Coagulation	(1W~40W, Load: 500 $\Omega$ )	It can be called Pinpoint Coagulation, it's relatively deep with the organization solidified layer, also the scope is comparatively narrow

### 3 Scope of Application

Applicable to Cutting and Coagulation of all kind of operations

### 4 Technical Parameter

Operating ambient temperature	5°C~40°C
Relative humidity range	80%
Pressure Range	86~106 kPa
Power Supply	220/110, 50/60 Hz
Operation Frequency	330~460 kHz

### 5 Transportation Data

Size: 21cm (L) \* 18cm (W) \* 9.5 (H) (Generator) ; 33 (L) \* 25 (w) \* 21cm (H) (Carton)  
Weight: 2kg (Net Weight) ; 3.5kg (Gross Weight)

### 6 Standard Accessories

1pcs of Power Cable, 2pcs of Disposable ESU Pencils, 2pcs Disposable ESU Pads,  
1pcs of Reusable cable for Disposable ESU Pad, 1pcs of Instruction Manual

