



# Shangyi Shadowless Operating Lamp SY02-LED3+3

## DESCRIPTION

1. Leading-edge lighting technology.
  2. Has great performance and outstanding energy efficiency,
  3. Illuminate demanding environments is becoming ever more common.
  4. Light emitting diodes have technological, economic and ecological advantages.
  5. Low heat generation, minimal power consumption and an almost unlimited life expectancy
- SY02-LED combines these characteristics. The result: Unparalleled performance.**

Type	Shangyi Shadowless Operating Lamp SY02-LED3+3
Cold light operation lamp Illuminance(Lux)	120000
Colour temperature(k)	4500±200
Colour reduction index(Ra)	96
Cold light operation lamp Illuminance depth(mm)	380/380
Total irradiance (W/m2)	200-350
Size of Cold light operation lamp field (mm)	200-350
Rated power of bulb(V/W)	3.2V/1W
Lamp head diameter(mm)	700/700
LED Service life of Cold light operation lamp illuminant	50000h
Power supply voltage	AC110V-240V(50HZ-60HZ)
Total power consumption(W)	120
Total LED bulb quantity(pcs)	222
Lowest height of Cold light operation lamp installation(mm)	2900

## FEATURES

1. **Efficient**  
Low power consumption and durable LEDs minimize the impact on our precious resources.
2. **Durable**  
50,000 LED life expectancy delivers tremendous product reliability and investment protection.
3. **Precise**  
Uniform, homogeneous lighting in surface and deep cavity procedures thanks to shadow control.
4. **Simple**  
Ergonomic, sterile operating concept allows simple and flexible controlling functions by surgical team.



2 piece in a carton: 104\*97\*44(cm)  
GW:36kg NW:27kg  
1 piece in a wooden case: 122\*42\*85(cm)  
GW:81kg NW:70kg

**AKL: 21603914144**